

# Boston BioProducts, Inc.

9 Short Street, Worcester, MA 01604

## Reagents for Biosciences

***Take 25% off your 1<sup>st</sup> order. Use code: 25BBP***

### Secondary Antibodies and Conjugates

Boston BioProducts secondary antibodies are prepared in goat using highly purified IgG as an immunogen. The IgG in the serum is immunoaffinity purified to remove all goat serum proteins including immunoglobulins. The affinity purified antibody is solid phased adsorbed with human IgG to ensure minimal cross reactivity. These secondary antibody conjugates are prepared and stored in preservatives.

<b>Product</b>	<b>Catalog #</b>	<b>Quantity</b>	<b>Price</b>
<b>Anti-mouse IgG (H+L, whole molecule)</b>	B-1105	0.5mg	\$45
	B-1110	1 mg	\$70
<b>Anti-mouse IgG (H+L)-Peroxidase</b> <i>Elisa 100,000 dilution</i> <i>Western Blot 1:10,000</i> <i>Immunostaining 1:200</i>	B-1115	0.5ml	\$60
	B-1120	1 ml	\$95
<b>Anti-mouse IgG (H+L)-FITC</b> Immunostaining 1:200	B-1125	0.5ml	\$60
	B-1130	1 ml	\$95
<b>Anti-mouse IgG (H+L)-Biotin</b> Elisa 100,000 dilution Immunostaining 1:200	B-1135	0.5ml	\$60
	B-1140	1 ml	\$95

<b>Product</b>	<b>Catalog #</b>	<b>Quantity</b>	<b>Price</b>
<b>Anti-rabbit IgG (H+L, whole molecule)</b>	B-1205	0.5mg	\$45
	B-1206	1 mg	\$70
<b>Anti-rabbit IgG (H+L)-Peroxidase</b> <i>ELISA 1:200,000 dilution</i> <i>Western Blot 1:10,000</i> <i>Immunostaining 1:200</i>	B-1215	0.5ml	\$60
	B-1220	1 ml	\$95
<b>Anti-rabbit IgG (H+L)-FITC</b> Immunostaining 1:200	B-1225	0.5ml	\$60
	B-1230	1 ml	\$95
<b>Anti-rabbit IgG (H+L)-Biotin</b> <i>ELISA 100.000 dilution</i> Immunostaining 1:200	B-1235	0.5ml	\$60
	B-1240	1 ml	\$95

***Take 25% off your 1<sup>st</sup> order. Use code: 25BBP***

Offer not valid if prohibited by institution