



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

## **SAFETY DATA SHEET**

### **SECTION 1 - PRODUCT AND COMPANY IDENTIFICATION**

- 1.1** PRODUCT NAME: 30% PEG 8000 Solution  
CATALOG #: BM-507
- 1.2** RECOMMENDED USE: Research use or further manufacturing
- 1.3** COMPANY IDENTIFICATION: Boston BioProducts, Inc.  
425 Fortune Blvd. Milford, MA 01757
- 1.4** EMERGENCY CONTACT: 1-800-535-5053  
TOLL FREE NUMBER: 1-888-881-8691

### **SECTION 2 – HAZARD(S) IDENTIFICATION**

#### **2.1 Classification of the substance or mixture**

**GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

Not hazardous

**Signal word:** None

**Hazard statement(s):** Not hazardous

#### **Precautionary Statements:**

|      |   |
|------|---|
| P101 | If medical advice is needed, have product container or label at hand. |
| p102 | Keep out of reach of children.  |
| P501 | Dispose of contents/ container to an approved waste disposal plant.   |

### **SECTION 3- COMPOSITION/INFORMATION ON INGREDIENTS**

**Composition:** 30 % (w/v) PEG 8000

| <b>Component</b> | <b>CAS Number</b> | <b>Concentration</b> |
|------------------|-------------------|----------------------|
| PEG 8000         | 25322-68-3        | 30 % (w/v)           |
| Water            | 7732-18-5         | Rest                 |

### **SECTION 4 - FIRST-AID MEASURES**

#### **4.1 Description of first aid measures**

**General advice:** Consult a doctor and show this safety data sheet.

- If inhaled:** Remove to fresh air and monitor breathing. If breathing becomes difficult, give oxygen. If breathing stops, give artificial respiration. Consult a doctor.
- In case of skin contact:** Immediately wash skin with copious amounts of soap and water for at least 15 minutes. Remove contaminated clothing and shoes and wash before reuse. Consult a doctor.
- In case of eye contact:** Flush with copious amounts of water for at least 15 minutes. Consult a doctor.



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

iv. **If swallowed:** Rinse mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Consult a doctor.

**4.2 Most important symptoms and effects, both acute and delayed:** To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

**4.3 Indication of immediate medical attention and special treatment needed:** Show this safety data sheet to the doctor in attendance. Immediate medical attention is required.

## **SECTION 5 - FIRE FIGHTING MEASURES**

### **5.1 Extinguishing media**

**Suitable extinguishing media:** Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide

**5.2 Special hazards arising from the substance or mixture:** Carbon oxides, Nitrogen oxides (NOx).

**5.3 Precautions for fire-fighters:** Wear suitable protective clothing to prevent contact with skin and eyes and self-contained breathing apparatus.

**5.4 Further information:** No data available

## **SECTION 6 - ACCIDENTAL RELEASE MEASURES**

**6.1 Personal precautions, protective equipment and emergency procedures:** Do not take action without suitable protective clothing - see section 8 of SDS. Evacuate personnel to safe areas. Ensure adequate ventilation. Avoid breathing vapors, mist, dust or gas.

**6.2 Environmental precautions:** Do not let product enter drains.

**6.3 Methods and materials for containment and cleaning up:** Cover spillage with suitable absorbent material. Using non-spark tools, sweep up material and place in an appropriate container. Decontaminate spill site with 10% caustic solution and ventilate area until after disposal is complete. Hold all material for appropriate disposal as described under section 13 of SDS.

**6.4 Reference to other sections:** For required PPE see section 8. For disposal see section 13.

## **SECTION 7 - HANDLING AND STORAGE**

**7.1 Precautions for safe handling:** Use in a chemical fume hood, with air supplied by an independent system. Avoid inhalation, contact with eyes, skin and clothing. Avoid the formation of dust and aerosols. Use in a well-ventilated area. Keep away from sources of ignition. Avoid prolonged or repeated exposure.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in cool, well-ventilated area. Keep away from direct sunlight. Keep container tightly sealed until ready for use. Keep in a dry place.

**7.3 Specific end use(s):** Use in a laboratory fume hood where possible.

## **SECTION 8 - EXPOSURE CONTROLS/PERSONAL PROTECTION**

### **8.1 Control parameters**

**Components with workplace control parameters:**

### **8.2 Exposure Controls**

- i. **Appropriate engineering controls:** Use in a fume hood where applicable. Ensure all engineering measures described under section 7 of SDS are in place. Ensure laboratory is equipped with a safety shower and eye wash station.

### **8.3 Personal protective equipment**



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

- i. **Eye/face protection:** Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).
- ii. **Skin protection:** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166
- iii. **Body Protection:** Wear appropriate protective clothing. Complete suit protecting against chemicals, the type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
- iv. **Respiratory Protection:** Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.
- v. **Control of environmental exposure:** Do not let product enter drains

## **SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES**

### **Information on basic physical and chemical properties**

|   |                           |                                  |                     |
|---|---------------------------|----------------------------------|---------------------|
| <b>Appearance</b>                                     | Clear, Colorless, Liquid  | <b>Vapor Pressure</b>            | < 0.01 hPa at 20 °C |
| <b>Odor</b>   | No Data Available         | <b>Vapor Density</b>             | No Data Available   |
| <b>Odor Threshold</b>                                 | No Data Available         | <b>Relative Density</b>          | 1.27 g/mL at 25 °C  |
| <b>pH</b>   | No Data Available         | <b>Water Solubility</b>          | No Data Available   |
| <b>Melting / Freezing Point</b>                       | 60 - 63 °C (140 - 145 °F) | <b>Partition Coefficient</b>     | No Data Available   |
| <b>Initial Boiling Point Range</b>                    | > 325 °C > 617 °F         | <b>Auto-Ignition Temperature</b> | 305 °C (581 °F)     |
| <b>Flash Point</b>                                    | No Data Available         | <b>Decomposition Temperature</b> | No Data Available   |
| <b>Evaporation Rate</b>                               | No Data Available         | <b>Viscosity</b>                 | No Data Available   |
| <b>Flammability (Solid, Gas)</b>                      | No Data Available         | <b>Explosive Properties</b>      | No Data Available   |
| <b>Upper / Lower Flammability or Explosive Limits</b> | No Data Available         | <b>Oxidizing Properties</b>      | No Data Available   |

## **SECTION 10 - STABILITY AND REACTIVITY**

**Stability:** Stable under recommended storage conditions

**Hazardous Decomposition Products/Hazardous Polymerization:** No data available

**Incompatibilities:** No data available

## **SECTION 11 - TOXICOLOGICAL INFORMATION**

### **Acute Toxicity**

| <b>Component</b> | <b>LD50 Oral</b>                 | <b>LD50 Dermal</b> | <b>LC50 Inhalation</b> |
|------------------|----------------------------------|--------------------|------------------------|
| PEG 8000         | LD50 Oral - Rat - > 50,000 mg/kg | -                  | -                      |



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

**Skin corrosion:** No data available

**Serious eye damage:** No data available

**Respiratory or skin sensitization:** No data available

**Carcinogenicity:** There are no known carcinogenic chemicals in this product

**Reproductive toxicity:** No data available

**Specific target organ toxicity – single exposure:** No data available

**Specific target organ toxicity – repeated exposure:** No data available

**Aspiration Hazard:** No data available

## **SECTION 12 - ECOLOGICAL INFORMATION**

**Eco Toxicity:** No data available

**Persistence and degradability/ Bioaccumulative potential/ Mobility in soil:** No data available

**Results of PBT and vPvB assessment:** No data available

**Endocrine disrupting properties:** No data available

**Other adverse effects:** No data available

**Chronic toxicity:** There are no known carcinogenic chemicals in this product

## **SECTION 13 - DISPOSAL CONSIDERATIONS**

Whatever cannot be saved for recovery or recycling should be managed in an appropriate and approved waste disposal facility. Processing, use or contamination of this product may change the waste management options. State and local disposal regulations may differ from federal disposal regulations. Dispose of container and unused contents in accordance with federal, state and local requirements.

## **SECTION 14 - TRANSPORT INFORMATION**

| <b>DOT</b>    | <b>TDG</b>    | <b>IATA</b>   | <b>IMDG/IMO</b> |
|---------------|---------------|---------------|-----------------|
| Not Regulated | Not Regulated | Not Regulated | Not Regulated   |

## **SECTION 15 - REGULATORY INFORMATION**

### **INTERNATIONAL INVENTORIES**

| <b>TSCA</b> | <b>DSL</b> | <b>NDSL</b> | <b>ELINCS</b> | <b>NLP</b> | <b>PICCS</b> | <b>ENCS</b> | <b>AICS</b> | <b>CHINA</b> | <b>KECL</b> |
|-------------|------------|-------------|---------------|------------|--------------|-------------|-------------|--------------|-------------|
| -           | -          | -           | -             | -          | -            | -           | -           | -            | -           |

### **USA FEDERAL REGULATIONS**

**SARA 313/SARA 311/312 HAZARDOUS CATEGORIZATION:** Not Hazard Class

**Clean Water Act/ Clean Air Act:** No Data Available

**Clean Air Act:** Not Applicable

**OSHA /TSCA 12(b):** Not Applicable

**California Proposition 65:** This product contains no Proposition 65 chemicals.

### **STATE RIGHT TO KNOW**

| <b>Massachusetts</b> | <b>New Jersey</b> | <b>Pennsylvania</b> | <b>Illinois</b> | <b>Rhode Island</b> |
|----------------------|-------------------|---------------------|-----------------|---------------------|
|----------------------|-------------------|---------------------|-----------------|---------------------|



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

|   |   |        |   |   |
|---|---|--------|---|---|
| - | - | Listed | - | - |
|---|---|--------|---|---|

**HMIS Rating**

| Health hazard | Chronic Health Hazard | Flammability | Physical Hazard |
|---------------|-----------------------|--------------|-----------------|
| 0             | 0                     | 0            | 0               |

**NFPA Rating**

| Health hazard | Fire Hazard | Reactivity Hazard |
|---------------|-------------|-------------------|
| 0             | 0           | 0                 |

**US DEPARTMENT OF TRANSPORTATION**

| REPORTABLE QUANTITY (RQ) | DOT MARINE POLLUTANT | DOT SEVER MARINE POLLUTANT |
|--------------------------|----------------------|----------------------------|
| No                       | No                   | No                         |

**U.S. DEPARTMENT OF HOMELAND SECURITY:** This product does not contain any DHS chemicals.

**Canada:** This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.

**This SDS complies with the requirements of Regulation (EC).**

**SECTION 16 - OTHER INFORMATION**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Boston BioProducts, Inc shall not be held liable for any damage resulting from handling or from contact with the above product. This product is sold for laboratory research and development purposes use only.

**Revised: 2025**