



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

## SAFETY DATA SHEET

### **SECTION 1 - PRODUCT AND COMPANY IDENTIFICATION**

- 1.1** PRODUCT NAME: Copper Chloride (0.1 M)  
CATALOG #: MT-301
- 1.2** RECOMMENDED USE: Research use or further manufacturing
- 1.3** COMPANY IDENTIFICATION: Boston BioProducts, Inc.  
425 Fortune Blvd. Milford, MA 01757
- 1.4** EMERGENCY CONTACT: 1-800-535-5053  
TOLL FREE NUMBER: 1-888-881-8691

### **SECTION 2 – HAZARD(S) IDENTIFICATION**

#### **2.1 Classification of the substance or mixture**

##### **GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

Acute toxicity, oral	Category 4
Acute toxicity, dermal	Category 4
Skin corrosion/ irritation	Category 2
Serious eye damage/eye irritation	Category 1
Hazardous to the aquatic environment, acute hazard	Category 1

#### **2.2 GHS Label elements, including Hazard and Precautionary Statement(s)**

##### **Pictogram**



**Signal word: Danger**

##### **Hazard statement(s):**

H302 +H312	Harmful if swallowed or in contact with skin.
H315	Causes skin irritation.
H318	Causes serious eye damage.
H410	Very toxic to aquatic life with long lasting effects.

##### **Precautionary Statements:**

P264	Wash skin thoroughly after handling.
P273	Avoid release to the environment.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352	IF ON SKIN: wash with plenty of soap and water.



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P501	Dispose of contents/ container to an approved waste disposal plant.

### **Section 3- COMPOSITION /INFORMATION ON INGREDIENTS**

**Composition:** 100 mM Copper (II) Chloride.

Component	CAS Number	Concentration
Copper (II) Chloride	7447-39-4	1.34 %
Water	7732-18-5	Rest

### **SECTION 4 - FIRST-AID MEASURES**

#### **4.1 Description of first aid measures**

**General advice:** Consult a doctor and show this safety data sheet.

- i. **If inhaled:** Remove to fresh air and monitor breathing. If breathing becomes difficult, give oxygen. If breathing stops, give artificial respiration. Consult a doctor.
- ii. **In case of skin contact:** Immediately wash skin with copious amounts of soap and water for at least 15 minutes. Remove contaminated clothing and shoes and wash before reuse. Consult a doctor.
- iii. **In case of eye contact:** Flush with copious amounts of water for at least 15 minutes. Consult a doctor.
- iv. **If swallowed:** Rinse mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Consult a doctor.

**4.2 Most important symptoms and effects, both acute and delayed:** To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

**4.3 Indication of immediate medical attention and special treatment needed:** Show this safety data sheet to the doctor in attendance. Immediate medical attention is required.

### **SECTION 5 - FIRE FIGHTING MEASURES**

#### **5.1 Extinguishing media**

**Suitable extinguishing media:** Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide

**5.2 Special hazards arising from the substance or mixture:** Carbon oxides, Nitrogen oxides (NO<sub>x</sub>).

**5.3 Precautions for fire-fighters:** Wear suitable protective clothing to prevent contact with skin and eyes and self-contained breathing apparatus.

#### **5.4 Further information:**

Health	Flammability	Instability	Physical Hazards
3	0	1	N/A

### **SECTION 6 - ACCIDENTAL RELEASE MEASURES**

**6.1 Personal precautions, protective equipment and emergency procedures:** Do not take action without suitable protective clothing - see section 8 of SDS. Evacuate personnel to safe areas. Ensure adequate ventilation. Avoid breathing vapors, mist, dust or gas.



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

**6.2 Environmental precautions:** Do not let product enter drains.

**6.3 Methods and materials for containment and cleaning up:** Cover spillage with suitable absorbent material. Using non-spark tools, sweep up material and place in an appropriate container. Decontaminate spill site with 10% caustic solution and ventilate area until after disposal is complete. Hold all material for appropriate disposal as described under section 13 of SDS.

**6.4 Reference to other sections:** For required PPE see section 8. For disposal see section 13.

## **SECTION 7 - HANDLING AND STORAGE**

**7.1 Precautions for safe handling:** Use in a chemical fume hood, with air supplied by an independent system. Avoid inhalation, contact with eyes, skin and clothing. Avoid the formation of dust and aerosols. Use in a well-ventilated area. Keep away from sources of ignition. Avoid prolonged or repeated exposure.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in cool, well-ventilated area. Keep away from direct sunlight. Keep container tightly sealed until ready for use. Keep in a dry place.

**7.3 Specific end use(s):** Use in a laboratory fume hood where possible.

## **SECTION 8 - EXPOSURE CONTROLS/PERSONAL PROTECTION**

### **8.1 Control parameters**

**Components with workplace control parameters:** Contains no substances with occupational exposure limit values.

#### **Exposure Guidelines**

<b>Component</b>	<b>ACGIH TVL</b>	<b>OSHA PEL</b>	<b>NIOSH IDLH</b>	<b>MEXICO OEL (TWA)</b>
Copper (II) Chloride	TWA: 1 mg/m <sup>3</sup>	-	IDLH: 100 mg/m <sup>3</sup> TWA: 1 mg/m <sup>3</sup>	-

### **Legend**

**ACGIH** – American Conference of Governmental Industrial Hygienists

**OSHA** – Occupational Safety and Health Administration

**NIOSH IDLH:** NIOSH National Institute for Occupational Safety and Health

### **8.2 Exposure Controls**

- i. **Appropriate engineering controls:** Use in a fume hood where applicable. Ensure all engineering measures described under section 7 of SDS are in place. Ensure laboratory is equipped with a safety shower and eye wash station.

### **8.3 Personal protective equipment**

- i. **Eye/face protection:** Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).
- ii. **Skin protection:** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166
- iii. **Body Protection:** Wear appropriate protective clothing. Complete suit protecting against chemicals, the type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
- iv. **Respiratory Protection:** Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

**v. Control of environmental exposure:** Do not let product enter drains

## **SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES**

### **Information on basic physical and chemical properties**

<b>Appearance</b>	Greenish to Blueish Liquid	<b>Vapor Pressure</b>	No Data Available
<b>Odor</b>	No Data Available	<b>Vapor Density</b>	No Data Available
<b>Odor Threshold</b>	No Data Available	<b>Relative Density</b>	No Data Available
<b>pH</b>	No Data Available	<b>Water Solubility</b>	No Data Available
<b>Melting / Freezing Point</b>	No Data Available	<b>Partition Coefficient</b>	No Data Available
<b>Initial Boiling Point Range</b>	No Data Available	<b>Auto-Ignition Temperature</b>	No Data Available
<b>Flash Point</b>	No Data Available	<b>Decomposition Temperature</b>	No Data Available
<b>Evaporation Rate</b>	No Data Available	<b>Viscosity</b>	No Data Available
<b>Flammability (Solid, Gas)</b>	No Data Available	<b>Explosive Properties</b>	No Data Available
<b>Upper / Lower Flammability or Explosive Limits</b>	No Data Available	<b>Oxidizing Properties</b>	No Data Available

## **SECTION 10 - STABILITY AND REACTIVITY**

**Stability:** Stable under recommended storage conditions

**Hazardous Decomposition Products/Hazardous Polymerization:** No data available

**Incompatibilities:** No data available

## **SECTION 11 - TOXICOLOGICAL INFORMATION**

### **Acute Toxicity**

#### **Component Information**

<b>Component</b>	<b>LD50 Oral</b>	<b>LD50 Dermal</b>	<b>LC50 Inhalation</b>
Copper (II) Chloride	LD50 Oral - Rat - 584 mg/kg	-	-

**Toxicologically Synergistic Products:** The environmental impact of this product has not been fully investigated

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

**Irritation:** Irritating to eyes

**Carcinogenicity:** The table below indicated whether each agency has any listed ingredient as a carcinogen.

<b>Component</b>	<b>CAS No</b>	<b>IARC</b>	<b>NTP</b>	<b>ACGIH</b>	<b>OSHA</b>	<b>Mexico</b>
-	-	-	-	-	-	-



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

## **SECTION 12 - ECOLOGICAL INFORMATION**

**Eco Toxicity:** The environmental impact of this product has not been fully investigated

Component	Freshwater Algae	Freshwater Fish	Microtox	Water Flea
Copper (II) Chloride	EC50: 0.12 - 0.2 mg/L/96h	LC50: 0.120-0.130 mg/L/96h (Carp) LC50: 0.9 mg/L/96h (Bluegill sunfish) LC50: 0.08 mg/L/96h (Rainbow trout)	-	EC50: 0.04 mg/L/48h

**Persistence and degradability:** No data available

**Bioaccumulative potential:** No data available

**Mobility in soil:** No data available

**Results of BPT and vPvB assessment:** No data available

### **Acute Toxicity**

Component	Log Pow
-	-

**Results of PBT and vPvB assessment:** No data available

**Other adverse effects:** No data available

**Chronic Toxicity:** There are no known carcinogenic chemicals in this product

## **SECTION 13 - DISPOSAL CONSIDERATIONS**

Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste.

Chemical waste generators must also consult local, regional, and national hazardous waste regulations to ensure complete and accurate classification.

Component	RCRA – U Series Wastes	RCRA – P series Wastes
-	-	-

Whatever cannot be saved for recovery or recycling should be managed in an appropriate and approved waste disposal facility. Processing, use or contamination of this product may change the waste management options. State and local disposal regulations may differ from federal disposal regulations. Dispose of container and unused contents in accordance with federal, state and local requirements.

## **SECTION 14 - TRANSPORT INFORMATION**

DOT	TDG	IATA	IMDG/IMO
UN Number: 2802 Name: Copper (II) Chloride Hazard Class: 8 Packing Group: III	UN Number: 2802 Name: Copper (II) Chloride Hazard Class: 8 Packing Group: III	UN Number: 2802 Name: Copper (II) Chloride Hazard Class: 8 Packing Group: III	UN Number: 2802 Name: Copper (II) Chloride Hazard Class: 8 Packing Group: III



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

**SECTION 15 - REGULATORY INFORMATION**

**INTERNATIONAL INVENTORIES**

TCA	DSL	NDSL	ELINCS	NLP	PICCS	ENCS	AICS	CHINA	KECL
Listed	Listed	-	Listed	Listed	Listed	Listed	Listed	Listed	Listed

**USA FEDERAL REGULATIONS**

**SARA 311/312 HAZARDOUS CATEGORIZATION**

Acute Health Hazard	Chronic Health Hazard	Fire Hazard	Sudden Release of Pressure Hazard	Reactive Hazard
Yes	No	No	No	Yes

**Clean Water Act/ Clean Air Act:** No Data Available

**Clean Air Act:** Not Applicable

**California Proposition 65:** This product contains no Proposition 65 chemicals.

**CERCLA:** This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response and Liability Act (CERCLA) (40 CRF 302)

Component	Hazardous Substance RQs	CERCLA EHS RQs
Copper II Chloride	10 lbs	10 lbs

**STATE RIGHT TO KNOW**

Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
Listed	Listed	Listed	-	-

**HMIS Rating**

Health hazard	Chronic Health Hazard	Flammability	Physical Hazard
3	0	0	0

**NFPA Rating**

Health hazard	Fire Hazard	Reactivity Hazard
3	0	1

**US DEPARTMENT OF TRANSPORTATION**

REPORTABLE QUANTITY (RQ)	DOT MARINE POLLUTANT	DOT SEVER MARINE POLLUTANT
Yes	No	No

**U.S. DEPARTMENT OF HOMELAND SECURITY:** This product does not contain any DHS chemicals.

**Canada:** This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

**This SDS complies with the requirements of Regulation (EC).**

#### **SECTION 16 - OTHER INFORMATION**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Boston BioProducts, Inc shall not be held liable for any damage resulting from handling or from contact with the above product. This product is sold for laboratory research and development purposes use only.

**Revised: 2025**